



Sports Brace

Sports Brace—Empowering People

TLSO BRACE

HCPCS: L0456 / L0457

SIZE: Universal



FEATURES:

- Easy to use universal sizing system
- Adjustable pull handles that allow for optimal compression
- Breathable fabric allows air flow to keep patients cool and dry over long periods of time

INDICATIONS: Spinal Stenosis, Spondylolistheses, Spondylolysis, chronic back pain, spinal instability, post-op disk surgery, Facet Syndrome; degenerative, herniated and bulging discs and post-op lumbar fusion.

HOW TO APPLY: Lay the brace out flat with the exterior viewing facing up. Open flap and remove both side panels. The size indications are marked on each panel. You can cut the brace with scissors at your appropriate size or simply fold the excess fabric and tuck it under the back of the panel. Attached the side panel to both the left and right side on the sizing panels. Position the back support around your waist and make sure the back panel is placed over the center of your lower back. Place the left panel securely over your abdomen and place the right panel over the left panel. Now that the brace is secure over your body, pull each of the compression straps simultaneously until proper compression is achieved.

CARE INSTRUCTIONS: Remove the rigid components. The rigid components can be wiped clean with a damp cloth. Soft goods can be hand washed in warm water with a mild detergent. Triple rinse and air dry. Avoid twisting the material to remove water.

Visit us online at www.sportsbrace.net